


## ·次性使用呼吸回路

Disposable breathing circuit





- 采用医用级高分子材料制造,无毒,无刺激
- 符合国际标准的接头,与国内外各种型号麻醉、 呼吸机相互匹配,连接牢固,密封性好
- L型接头配有检测口,便于取样和检测
- 集水杯可收集储存管路内的冷凝水
- 单个病人一次性使用,避免交叉感染
- 规格:成人,儿童



-次性使用呼吸回路(波纹管) Corrugated type

- Made of medical grade polymer material, non-toxic and no stimulation
- Connector accord with ISO standard, match with all Anesthesia & Breathing machines, secure connection, reliably sealing
- L type connector with detection port for sampling and testing
- Water traps can collect condensate water from breath circuit.
- Single use and avoid cross infection
- Size: Adult and Pediatric

管长 / Tube length

1.2m

1.6m

1.8m

可依客户定制 Customized

## 零部件介绍 / Components introduction





## 双垂直接口,使积水杯随时保持垂直状态 最大化使用积水杯的功能,减少管内积水



The hardness is high, good anti-distortion not easy to bend, the inner wall is smooth without hydrops.





## 广东海鸥医疗器械股份有限公司 GUANGBONG HAIBU MEDICAL APPARATUS CO., LTD. TEL: +86-663-2758666 FAX: +86-663-2905999

广东省普宁市流沙北南园工业区 Liusha, Puning City, Guangdong, China







美国FDA认证